

# SCOOTERPAC FEATHER PROPEL

SA2020

\$1,590



**LIGHTWEIGHT  
COMPACT DESIGN**



**SOLID PUNCTURE  
PROOF TYRES**



**KG**

**113KG WEIGHT  
CAPACITY**



# SCOOTERPAC

## FEATHER PROPEL SA2020

Weighing just 7.2kg, the Feather Propel is one of the lightest wheelchairs on the planet. Taking advantage of the latest materials, the Feather Propel sets a new standard. Practically anyone can lift this wheelchair into their car plus it is super easy to propel and push. Unlike other wheelchairs, Propel offers optional extras such as elevating leg rests for superior comfort and support.

Transporting the Feather Propel is a breeze. The chair folds down into a compact package. Simply fold down the backrest and then pull up the middle of the seat.

### SPECIFICATIONS

<b>Adjustable Armrests</b>	Yes	<b>Seat Dimensions (wxd)</b>	457mm x 406mm
<b>Adjustable Backrest</b>	No	<b>Seat height from floor</b>	483mm
<b>Adjustable Seat Height</b>	No	<b>Seat Height from footplate</b>	TBC
<b>Adjustable Handle Height</b>	No	<b>Seat padding</b>	Medium
<b>Armrests</b>	Yes 229mm	<b>Tyres</b>	Solid, Polyurethane
<b>Dimensions (L x W x H)</b>	889 x 762 x 635	<b>Weight Capacity</b>	113kg
<b>Dimensions when folded</b>	381 x 737 x 711	<b>Total Weight</b>	7.2kg
<b>Distance between armrests</b>	TBC	<b>Backrest height</b>	406mm
<b>Drive Wheel</b>	559mm	<b>Main Material</b>	Aluminium
<b>Front Wheel</b>	152mm	<b>Brakes</b>	attendant handles & push to lock
<b>Heaviest part</b>	6.1kg	<b>Footplate options</b>	Swing-away
<b>Number of components</b>	5	<b>Shipping weight</b>	