

MIO AVANZO ROLLATOR

PRODUCT CODE: CA891L



WEIGHS JUST
6KG



SOLID PUNCTURE
PROOF WHEELS



LARGE STORAGE
BAG & BASKET



MIO AVANZO ROLLATOR

PRODUCT CODE: CA891L

Introducing the Avanzo Rollator with a Large Storage Basket, a versatile and reliable mobility aid designed to enhance your daily activities and provide you with a convenient and secure way to carry your belongings. The Avanzo Rollator combines functionality, durability, and style to offer you an exceptional walking experience.

Its sturdy frame is constructed from high-quality materials, ensuring long-lasting performance and stability. The rollator is engineered to support a wide range of individuals, providing a comfortable and secure walking assistance for those with mobility challenges or anyone seeking additional support while on the move.

SPECIFICATIONS

Weight Capacity	100kg	Total Weight	6kg
Total Length	630mm	Folded Dimensions	630mm x 250mm x 820mm
Total Width	590mm	Heaviest Part	5.1kg
Handle Height	820-980mm	Storage	Bag and Basket