

SOLAX OMNI WHEEL POWER CHAIR

FOLDING POWER CHAIR



**JOYSTICK
CONTROLLER**



**PUNCTURE PROOF
TYRES**



**LIGHTWEIGHT
COMPACT DESIGN**



SOLAX OMNI-WHEEL POWER CHAIR **FOLDING POWER CHAIR**

Discover the Solax Omni Wheel Automatic Folding Electric Wheelchair, a state-of-the-art mobility solution that effortlessly combines convenience and reliability. This wheelchair features an automatic folding mechanism, streamlining transport and storage for users on the move. With a robust design, it supports a weight capacity of up to 125kg, ensuring a secure and comfortable experience.

SPECIFICATIONS

Adjustable Armrests	No	Maximum Speed	7km/h
Adjustable Backrest	No	Motor	Brushless
Adjustable Seat Height	No	Number of components	1 + battery (x2)
Anti-Tip Mechanism	Yes	Seat Dimensions (w x d)	450 x 430mm
Battery	24V 2Ah Lithium	Seat height from floor	535mm
Battery Range	24km	Seat Height from footplate	430mm
Battery Weight	1.3kg	Seat padding	Medium
Dimensions	1030 x 600 x 875mm	Slope ability	8°
Dimensions when folded	825 x 600 x 410mm	Suspension	Yes
Distance between armrests	450mm	Swivel seat	No
Drive Wheel	260mm	Turning Radius	600mm
Front Wheel	260mm	Tyres	Puncture Proof Omni-wheels
Ground Clearance	105mm	Weight Capacity	125kg
Headrest	No	Weight with Battery	49.3
Heaviest part	48kg	Weight without battery	48kg
Lights	Yes	Backrest height	410mm
Lithium Battery	Yes	Available Sizes	Standard